

Evaluating the efficacy of the MicroO₂ mandibular repositioning device for the treatment of obstructive sleep apnea

Jerry Hu, DDS, David Kuhns, PhD, Sung Kim, BS, and Len Liptak, MBA

Obstructive Sleep Apnea (OSA) is a rapidly emerging public health concern that enables dentists to expand the scope of their services. This single-center study offers a first look at the clinical efficacy of the new MicroO₂ Sleep Apnea Device.

BACKGROUND

The American Academy of Dental Sleep Medicine notes that over 18 million Americans have sleep apnea.¹ The long-term effects of sleep loss and sleep disorders have been associated with increased risk of hypertension, diabetes, obesity, depression, heart attack, and stroke.

Medical and dental professionals are becoming increasingly aware of the complexity, difficulty, and seriousness of treating sleep issues, making the concurrent use of gold standard treatments (CPAP, BIPAP, APAP, ASV) and dental mandibular repositioning devices (MRD's) a viable and often prudent treatment approach.

STUDY DESIGN

The goal of this study was to evaluate the clinical efficacy of the new MicroO₂ Mandibular Repositioning Device (MicroDental Laboratories, Dublin, CA). For this study, clinical efficacy is defined according to three success thresholds that are commonly utilized in the field of dental sleep medicine.² The first success threshold is a post treatment AHI score of less than 5. The second is a post treatment AHI score of less than 10. The third is a post treatment AHI score improvement (reduction) of 50% or better, if not exhibiting an AHI score below 10.

This single-center study is comprised of either patients recently diagnosed with obstructive sleep apnea or existing patients who required a new MRD. Each patient was treated with the MicroO₂ Sleep Apnea Device.

STUDY RESULTS

For the six participants who completed the protocol at the time of authoring this report, 100% achieved at least one of the three success thresholds.

CONCLUSION

The findings of this single-center study of six participants confirm that the new MicroO₂ MRD is a viable option for treating patients who are diagnosed with obstructive sleep apnea. Further research is required to evaluate the potential

incremental clinical benefits of this new MRD with respect to patient compliance, treatment efficiency, and mitigation of MRD side effects.


ADDITIONAL INFORMATION

The Las Vegas Institute, The American Academy of Dental Sleep Medicine, and The Pankey Institute are amongst the leading providers of continuing education courses on the practice of dental sleep medicine. For more information, contact these institutions and organizations, or contact MicroDental Laboratories.

REFERENCES


- 1 www.AADSM.org
- 2 Ferguson KA; Cartwright R; Rogers R et al. Oral Appliances for Snoring and Obstructive Sleep Apnea: A Review. *SLEEP* 2006;29(2): 244-262.

Micro₂MRD Therapy Pre-Post AHI Score Analysis, by Participant




NEW

MicrO₂™ SLEEP DEVICE. MORE SLEEP. LESS HASSLE.



LINGUAL-FREE DESIGN
Offers more tongue space

DIGITAL PRECISION
Provides retention by mirroring dentition



TWIN 90° ANGLE MECHANISM-FREE POSTS
Keeps jaw forward at night

CONTOURED LIP AND CHEEK BORDERS
Provides comfort

The lingual-free MicrO₂ is engineered to treat your patients with greater comfort, predictability and simplicity.

What are recognized leaders in dental sleep medicine saying about their initial MicrO₂ cases?

"The CAD/CAM design optimizes outcome but minimizes bulk—a true break through advancement in the treatment of snoring and sleep apnea." **DR. ANNE-MAREE COLE, BDS, LVIM**

"We found that with the added tongue room the mandible did not have to be brought as far forward. The patient felt overall it was less bulky and more comfortable."

DR. NANCY ADDY, DDS, DIPLOMAT AADSM

"With its airway constriction preventing fins and unobtrusive size, this device is a winner."

DR. WILLIAM G. DICKERSON, DDS, FAACD, LVIM

"It's small, tough, and reliable. Having no tiny adjustment mechanism is helpful for some patients—the easy way to bring the jaw forward and keep the airway open." **DR. STEVE CARSTENSEN, DDS**

"The MicrO₂ sleep appliance is so simple, effective and easy for patients to wear."

DR. MARK T. MURPHY, DDS, FAGD

Discover if MicrO₂ is simply a great fit for your practice.

Contact an Authorized MicroDental Network Lab for more information and to pre-schedule a case.

| | |
|--|---------------------|
| USA, NATIONWIDE | 800.229.0936 |
| Chicago, IL (MicroDental Chicago)..... | 888.622.7835 |
| Cincinnati, OH (Mitch Lab)..... | 800.767.5303 |
| Dallas, TX (Westbrook Lab)..... | 800.718.3384 |
| Detroit, MI (Ward Lab)..... | 800.833.3865 |
| Hermiston, OR (Aesthetic Designs Lab)..... | 541.567.1667 |
| Indianapolis, IN (Mitch Lab)..... | 800.767.5303 |
| Lansing, MI (Precision Lab)..... | 800.292.0855 |
| Las Vegas, NV (MicroDental at LVI)..... | 800.933.6838 |
| Louisville, KY (ADLab)..... | 800.456.1292 |
| New York, NY (TetraDynamics)..... | 800.877.8271 |

| | |
|--|---------------------|
| Raleigh, NC (MicroDental Research Triangle)..... | 800.840.2651 |
| Rochester, NY (CQCLab)..... | 800.724.1058 |
| Seattle, WA (DiMartino Lab)..... | 800.562.0300 |
| Saint Paul, MN (Hermanson Lemke Lab)..... | 800.328.9648 |
| Salt Lake City, UT (MicroDental Salt Lake)..... | 888.344.9991 |
| San Francisco Area, CA (MicroDental)..... | 800.229.0936 |
| St. Louis, MO (Donnell Lab)..... | 800.344.7866 |
| CANADA, NATIONWIDE | 866.222.0035 |
| Ottawa, ON (MicroDental Ottawa)..... | 866.222.0035 |
| Hamilton, ON (MicroDental Hamilton)..... | 800.263.7250 |
| Vancouver, BC (MicroDental Vancouver)..... | 800.561.0926 |

MicrO₂
SLEEP AND SNORE DEVICE

FDA | PATENT PENDING DESIGN
CLEARED



SMILES MATTER™
MicroDental.com

MICRODENTAL
LABORATORIES